Navodila za uporabo – ActiMaris® NASAL pršilo

1. Opis proizvoda

ActiMaris® NASAL pršilo je hipertonična, alkalno ionizirana raztopina morske soli (1,2% NaCl), natrijevega hipoklorita NaOCl (0,04%) in HOCl hipoklorove kisline (0,004%) in deluje na bio-fizični osnovi. Morska sol in natrijev hipoklorit ustvarjata protivnetni in protimikrobni učinek.

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Medizinprodukt: Klasse I, CE

EC REP

Adriatic Health Factory d.o.o.
Nikole Tesle 17
44000 Sisak
CROATIA

 PHARM

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Z-Pharm d.o.o.
Dunajska cesta 181
1000 Ljubljana

CE

Instructions for use – ActiMaris® NASAL Spray

1. Product profile

ActiMaris® NASAL Spray is a hypertonic alkaline ionised seawater solution with sea salt (1.2 %), sodium hypochlorite NaOCl (0.04 %) and hypochlorous acid HOCl (0.004 %) and acts on a bio-physical basis. The sea salt and sodium hypochlorite create an anti-inflammatory and antimicrobial effect.

ActiMaris® NASAL Spray moisturises and cleanses in interaction with the mucous membrane. ActiMaris® NASAL Spray also supports the physiological milieu on the surface of the mucous membrane through its composition of minerals, salts and trace elements. Due to its hypertonicity/hypermolarity, the solution has a decongestant effect, dissolves plaque and nasal secretions and frees nasal breathing. Thanks to its cleansing effect and dissolving of secretions, ActiMaris® NASAL Spray improves well-being and alleviates the typical symptoms of the above-mentioned indications. ActiMaris® NASAL Spray is both cell-friendly and antimicrobially efficient and fast. The sea salt and sodium hypochlorite are natural preservatives and guarantee a long stability of the ActiMaris® NASAL Spray solution.

ActiMaris® NASAL Spray contains:

- no cortisone
- no xylometazoline
- no alcohol
- no antibiotic

2. Fields of application/indications

ActiMaris® NASAL Spray is suitable as an additional and supporting therapy of:

- Rhinitis/cold/ coryza
- Influenza infections
- Acute and chronic inflammation of the nasal mucous membrane (rhinitis acuta et chronica)
- Acute and chronic sinusitis (rhinosinusitis acuta et chronica)
- Rhinitis medicamentosa after abuse of sympathomimetic nasal sprays
- Allergic rhinitis (also prophylactic)
- Tube catarrh/middle ear infections
- Post-operative (after nose and sinus surgeries)
- Prophylaxis of multi-resistant germs MRSA/VRE
- Prophylaxis of inflammation in the lower respiratory tract (reduces post-nasal drip)

3. Application instructions and dosage

If the head is bent slightly forward, the spray reaches the middle and lower nasal passage. In this way the covered surface area of the mucous membrane is increased and a better effect is achieved.

Adults and children from 2 years:

Recommended dosage of ActiMaris® NASAL Spray: 2 sprays in each nostril 4–5 times daily.

Can be used during pregnancy and lactation.

Toddlers (under 2 years) and infants:

1–2 sprays in each nostril 2–3 times daily as needed.

Infants and toddlers should not be placed on their backs, but instead on the side when ActiMaris® NASAL Spray is applied. The head should be lifted slightly. Always consult your paediatrician before using the product for infants and toddlers.

For more severe endonasal infections and mucous membrane coatings, ActiMaris® NASAL Spray can be combined with ActiMaris® wound gel 2–3 times daily. At the recommended dosage, the content of 20 ml (100µl per spray burst) is sufficient for a therapy of between 7 and 10 days (approx. 200 spray bursts). When using ActiMaris® NASAL Spray, an overdose is not possible, nor is there any systemic effect. The nasal solution disintegrates into water, salt and oxygen after contact and reaction with the mucous membrane. In the same way a dependence on/habituation to the nasal solution is not possible.

4. Tissue tolerance and biocompatibility

ActiMaris® NASAL Spray is suitable for long-term (>24h) and repeated use in irritated mucous membrane processes, as tissue tolerance and biocompatibility have been proven also with repeated or multiple applications. The application is almost painless and is generally well tolerated by allergy sufferers.

5. Side effects and contraindications

So far, no side effects are known.

Contraindication in case of allergy to ingredients (not observed so far).

6. Interactions

When used correctly, no interactions are to be expected – even with necessary systemic antibiotic therapies. A local combination therapy with another endonasal preparation is not necessary or recommended due to its optimal effect in the stated indications. Please consult your doctor for safety reasons, if you want to systematically use other drugs in combination.

7. General safety instructions

Keep out of the reach of children! Please avoid direct contact with coloured fabrics, as the contained sodium hypochlorite can have an oxidising or bleaching effect! Rinse immediately with water after contact with coloured fabrics.

8. Ingredients and technical data

Water, sea salt, sodium hypochlorite.

Clear, colourless, alkaline and hypertonic sea salt solution with characteristic odour; contains 1.2 % sea salt (Sal Maris), 0.04 % NaOCl and 0.004% HOCl.

9. Storage instructions, shelf life, period of use

ActiMaris® NASAL Spray should be stored at room temperature between +5 °C and +25 °C. Protect the spray from direct sunlight and/or heat.

The shelf life is 30 months. The current expiration date is printed on the packaging.

After opening the spray, the period of use is 3 months.

ActiMaris AG

natural medical solutions

Manufacturer:

ActiMaris AG
Sandgrube 29
9050 Appenzell (SWITZERLAND)
info@actimaris.com
Phone +41 71 505 75 25

 Batch designation/batch number
 Order number
 Attention, observe instructions for use
 Can be used until
 Storage instructions (temperature)
 Protect from direct sunlight
 Manufacturer

Information status: 10/2020 – V1.2

Medicinal product: Class I, CE

 PHARM

Distributor Slovenia:

Z-Pharm d.o.o.
Dunajska cesta 181
1000 Ljubljana

EC REP

Adriatic Health Factory d.o.o.
Nikole Tesle 17
44000 Sisak
CROATIA

CE