

CONTACT PLUS® NÉO Toilet raiser



NEW DESIGN

XXL toilet



2-YEAR WARRANTY
MAX. USER'S WEIGHT: 185 KG

The new generation of **Contact plus Néo** toilet raisers brings a comfort never equalled on the market thanks to the Herdegen Softinject technology. With its double fixing lips, **Contact Plus Néo** quickly and securely fits on all standard toilet bowls to raise their height by 11 cm.


Materials: injected HR foam, PVC skin
Cleaning: non-abrasive cleaner, dry with a soft cloth
Autoclavable (5min at 90°C and 80% humidity)

- > Tested and certified up to 185 kg
- > High-resilience foam and soft PVC skin
- > Global innovation: unmatched comfort and safety



Video presentation

Packaging:
Pieces per box: 1
Pallet: 66



Unmatched comfort and safety



Easy to install without screws or tools

TECHNICAL FEATURES

	Contact plus ® Néo
Overall width	37 cm
Seat opening	21*26 cm
Toilet elevation	11 cm
Overall height	16 cm
Overall depth	41 cm
Weight	2 kg



Retail box

IDENTIFICATION CODES

Code	Item	EAN Code	HS Code	GMDN	CND
500200n	Contact plus ® Néo	3701127704426	392490009	091218	Y091215

