

# CLIPPER® IV & V Toilet raisers



Clipper® V

The **Clipper® IV & V** raise the toilet seat by 11 cm and can be installed without any tools. Their removable armrests are padded and enable the patient to stand up and sit down easier. **Clipper V** (500421) is equipped with a lid for optimal hygiene.


Materials: 100% polypropylene  
 Cleaning: non-abrasive cleaner, dry with a soft cloth, autoclavable

- > Tested and certified up to 100 kg
- > Smooth and shiny finish
- > Easy fixation without any tools or screws
- > Removable armrests
- > Maximum comfort and security

Video presentation



Packaging:  
 Pieces per box: 1  
 Pallet: 60



## TECHNICAL FEATURES

	Clipper IV®	Clipper V®
Overall width	66 cm	
Seat opening	21*25 cm	
Toilet elevation	11 cm	
Overall depth	42 cm	
Weight	2.7 kg	3.1 kg



Retail box

## IDENTIFICATION CODES

Code	Item	EAN Code	HS Code	GMDN	CND
500420	Clipper IV®	3701127709056			
500421	Clipper V® with lid	3701127709063	392490009	091218	Y091212

