

BATHTUB HANDLE «FIXY»



- ✓ Locking system easy to use
- ✓ Treated steel frame, epoxy painting
- ✓ Stabilization screw
- ✓ Installation without tools



DESCRIPTION

This bathtub handle is to be placed by tightening the screw on the bathtub's side (Thickness min 10 and 18,5 cm max) The handle has a useful height of 38cm to allow the user to span the bathtub using the bathtub's point of support. Please make sure to lock properly the handle before use.

Max user's weight 100 kg

Warranty : 2 years

Ref	Designation	Total height	Handle height	Total width	Handle width	Depth	Weight
550805	Bathtub handle Fixy	50 cm	38 cm	15,2 cm	14,5 cm	19,5 cm	2,96 kg